

Transfer Board

HCPCS:



- For ease of transferring to and from a wheelchair
- Made of sturdy Baltic birch with a clear lacquer finish
- Ends tapered, corners rounded
- Weight Capacity 440 lbs

| PRODUCTS | Description |
|----------|---------------------------------------|
| RTL6045 | 30" Transfer Board w/Handles,1EA,2/CS |

| SPECIFICATIONS | Description |
|-----------------|------------------|
| Warranty | Limited 6 Months |
| Weight Capacity | 440 lb |