

## Bariatric Aluminum Two Button Folding Walker with Wheels

HCPCS:




FIGURE A



- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design when compared to standard walkers
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability
- Vinyl-contoured hand grip

| PRODUCTS               | Description                       |
|------------------------|-----------------------------------|
| 10220-1WW              | Walker, AdltHD, 2Btn w/Whls, 1/cs |
| SPECIFICATIONS         | Description                       |
| Warranty               | Limited Lifetime                  |
| Weight Capacity        | 500 lb                            |
| Height                 | 32" - 39"                         |
| Carton Shipping Weight | 24 lb                             |
| Overall Width          | 27"                               |
| Weight                 | 11 lb                             |