

## Steerable Knee Walker

HCPCS: E0118



- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- New, two-piece padded leg pad provides comfort
- Deluxe braking system
- Removable front basket
- Durable, steel frame with silver vein finish
- Tool free height adjustment
- 7.5" casters are ideal for indoor/outdoor use

| PRODUCTS               | Description                                        |
|------------------------|----------------------------------------------------|
| 796                    | Steerable Knee Walker,Split KneePad,1/EA           |
| SPECIFICATIONS         | Description                                        |
| Warranty               | Limited Lifetime, Lifetime Warranty on Brake Cable |
| Handle Height          | 31.5"-40.25"                                       |
| Weight Capacity        | 350 lb                                             |
| Carton Shipping Weight | 27 lb                                              |
| Knee Pad Dimensions    | 7" (W) x 15" (D)                                   |
| Knee Pad Height        | 16.75"-20.5"                                       |
| Weight                 | 21 lb                                              |
| Wheel Size             | 7.5"                                               |