

Steerable Knee Walker

HCPCS: E0118



- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- New, two-piece padded leg pad provides comfort
- Deluxe braking system
- Removable front basket
- Durable, steel frame with silver vein finish
- Tool free height adjustment
- 7.5" casters are ideal for indoor/outdoor use

PRODUCTS	Description
796	Steerable Knee Walker,Split KneePad,1/EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime, Lifetime Warranty on Brake Cable
Handle Height	31.5"-40.25"
Weight Capacity	350 lb
Carton Shipping Weight	27 lb
Knee Pad Dimensions	7" (W) x 15" (D)
Knee Pad Height	16.75"-20.5"
Weight	21 lb
Wheel Size	7.5"