

Viper Plus Reclining Wheelchair

HCPCS: K0004, E1226-1,K0004, E1226-1



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Patented built-in seat rail extension and extendable upholstery easily adjusts seat depth from 16" or 18"
- Attendant wheel locks in a convenient and easily accessible location
- Spreader bar acts as a back stabilizer and a convenient push handle for the attendant
- Offset rear wheel position adds additional stability and equally distributes the patient's weight while in the various reclining positions
- Hydraulic reclining mechanism allows for adjustments up to 180°
- Durable, flame-retardant nylon upholstery resists mildew and bacteria
- Carbon steel frame with double-coated chrome for an attractive, chip-proof, maintainable finish
- Composite Mag-style wheels are lightweight and maintenance free also features 8" front casters
- Nylon, padded armrests provide added user comfort, and feature the flip-back style for easy transfers
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Standard dual axle provides easy adjustment of seat-to-floor height to hemi-level
- Front caster forks are adjustable in three positions
- Comes standard with swing-away elevating leg rests
- Composite foot plates are durable and lightweight
- Polyurethane rear tires and casters offer superior performance, smoother ride and are lightweight

PRODUCTS

Description

PLA416RBDDA	ViperGT,16"Recl,FlpBackRemDsk,1/ea
PLA416RBDFA	ViperGT,16"Recl,FlpBackRemFull,1/ea
PLA418RBDDA	ViperGT,18"Recl,FlpBackRemDesk,1/ea
PLA418RBDFA	ViperGT,18"Recl,FlpBackRemFull,1/ea
PLA420RBDDA	ViperGT,20"Recl,FlpBackRemDesk,1/ea
PLA420RBDFA	ViperGT,20"Recl,FlpBackRemFull,1/ea

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Product Weight	47 lb 49 lb 51 lb
Product Dimensions	42" (L) w/ riggings x 24" 26.5" 28" (W) x 50" (H)
Seat Dimensions	16" 18" 20" (W) x 16" - 18" (D)
Weight Capacity	300 lb
Seat To Floor Height	17.5" 19.5"
Carton Shipping Weight	53 lb 55 lb 57 lb
Folded Dimensions	42" (L) w/ riggings x 12.5"(W) x 50" (H)
Back Height	23" 33" w/ headrest

ACCESSORIES	Description
STDS6008-1	Carry Bag,W/C,Oxygen,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS834	OverheadAnti- Theft Bar,Single,1/ea, EA
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS831	Heel Loops, 1pr, BX
WASL	Limb Support,Left,1/cs, CS
STDS801	Brake Handle Ext.,8", 1/ea, EA
WASR	Limb Support, Right,1/cs, CS
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, EA
PARTS	Description
STDS4S4803	HeadrestAssembly, 18", CS Recliner, 1/, EA
STDS4S4827	Back Upholstery, 18", CS Recliner, 1/e, EA
STDS4S4802	HeadrestAssembly, 16", CS Recliner, 1/, EA
STDS4S4801	Headrest Pillow, CS Recliner, 1/ea, EA
STDS4S4808	Backrest Cane, CS Recliner, 1/ea, EA
STDS4S4826	Back Upholstery, 16", CS Recliner, 1/e, EA
STDS4S4828	Back Upholstery, 20", CS Recliner, 1/e, EA
TRHG	Hand Grips/Universal,1pr/bx, BX
STDS4Y4702L	ReleaseHndle w/Cable,SSP18RBDDA,LFT,1/ea, EA
STDS4Y4702R	ReleaseHndle w/Cable,SSP18RBDDA,RHT,1/ea, EA
STDS4Y4715	Anti Folding ADJ Knob,SSP18RBDDA,1/ea, EA
STDS4S4809	SpreaderBar,16",CS,Sentra,Viper Rec,1/ea, EA
STDS4S4810	SpreaderBar,18",CS,Sentra,ViperRec,1/ea, EA
STDS4S4811	SpreaderBar,20",CS,Sentra,Viper,Rec,1/ea, EA
STDS3006	Cable,Back Release,Reclining w/c, 1/ea, EA
STDS4Y4704L	Piston,allCompnts,SSPRBDDA,Left,1/ea, EA
STDS4Y4704R	Piston,allCompnts,SSPRBDDA,Right,1/ea, EA
STDS3J2400	Handrim,Cruiser3,1/ea, EA
STDS004609-9	Hand Rim Screw Set of 9,1Bag, BAG
STDS4S4830	RearWheelAxleBolt/Nut/Spqr,CSReclnr,1ea, EA