

Blue Streak Wheelchair

HCPCS: K0001-1,K0001-1



- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- Solid Polyurethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap (-SF models only)
- Comes with swing-away footrests or elevating leg rests

PRODUCTS	Description
BLS16FBD-ELR	BlueStreak16"D/L FlipBackArms,Elev,1/ea
BLS16FBD-SF	BlueStreak16"D/L FlipBackArms,SwgFt,1/ea
BLS18FBD-ELR	BlueStreak18"D/L FlipBackArms,Elev,1/cs
BLS18FBD-SF	BlueStreak18"D/L FlipBackArms,SwgFt,1/cs
BLS20FBD-ELR	BlueStreak 20"D/FlipBackArms,EivLeg,1/ea
BLS20FBD-SF	BlueStreak 20"D/FlipBackArms,SwgFt,1/ea

SPECIFICATIONS	Description
Warranty	Limited 3 Year
Product Weight	41.25 lb
Weight Capacity	300 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	47 lb
Back Height	16"

ACCESSORIES	Description
STDS833N	Rear Anti Tipper w.wheels SSP&CS,1/pr, PR
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS831	Heel Loops, 1pr, BX