

Airgo Comfort-Plus XC Premium Transport Chair

HCPCS:



- Designed to be as pleasurable as it is practical
- Angled handlebars provide more room for one's back, while caregivers will appreciate the comfortable ergo-form handgrips with convenient lockable hand brakes for increased safety
- Unique footrests can be conveniently locked to the side of the frame for easy storage
- Detachable flip back armrests also make it easier to pull up to a table or desk
- Heel-loops and seat belt included as standard features for increased safety
- Large 12" rear wheels with premium bearings for smoother rolling and increased comfort
- Ergo-form handgrips comfortably fit your palm and provide natural brake access
- Handgrip mounted brake levers for increased control
- Anatomically correct seat-to-armrest height for a naturally erect posture
- Detachable full flip-back armrests for effortless transfers
- Fold down back for more compact transport and storage
- Convenient backrest pouch for discreet storage

PRODUCTS

Description

700-830	Airgo Comfort-Plus XC TC, Black, 1 EA
---------	---------------------------------------

SPECIFICATIONS

Description

Warranty	Limited lifetime on frame 90 days on wearable components
Product Weight	21 lbs
Weight Capacity	300 lbs
Seat To Floor Height	19"
Carton Shipping Weight	23 lbs