

## 1 Drawer Bedside Cabinet

HCPCS:



- High Pressure wood grained Laminate Cabinet finish, Buried wood edges with High Pressure Laminate Top, Drawer front and door feature raised real wood detail, Solid wood feet
- Available in six most popular wood grain finishes
- Drawer slides are epoxy coated and have 100 lbs rating (Figure A)
- European style hinges provide 110° accessible opening and self closing feature (Figure B)
- Ventilated back allows air flow throughout cabinet interior (Figure C)
- Warranty: Seven years

### PRODUCTS

RS-BDSC-11

### Description

Riverside Bedside 1dr/1drwr