

3 Drawer Bedside Cabinet

HCPCS:

- High Pressure wood grained Laminate Cabinet finish, Buried wood edges with High Pressure Laminate Top, Drawer front and door feature raised real wood detail, Solid wood feet
- Available in six most popular wood grain finishes
- Drawer slides are epoxy coated and have 100 lbs rating (Figure A)
- Ventilated back allows air flow throughout cabinet interior (Figure B)
- Available in Arched and Flat Toe Kicks (Figure C)

PRODUCTS

Description

RS-BDSC-03

Riverside Bedside 3drwr
