

## 4 Drawer Chest

HCPCS:

- Top & Sides feature High Pressure Thermo fused Melamine with matching vinyl edge on top
- Available in nine most popular wood grain finishes
- Drawer slides are epoxy coated and have 100 lbs rating (Figure A)
- Ventilated back allows air flow throughout cabinet interior (Figure B)
- Available in Arched and Flat Toe Kicks (Figure C)
- Warranty: Seven years

<b>PRODUCTS</b>	<b>Description</b>
BW-CHST-04	Brownlee Woods 4 Drawer Chest