

## Oslo 3 Drawer Bedside Cabinet

HCPCS:



- Top High Pressure Thermo Fused Laminate
- Full Extension, Ball Bearing Drawer slides with rubber stops provide smooth operation, quiet closure
- Locking top drawer is a standard
- Contemporary drawer handles with brushed nickel finish
- Ventilated back allows air flow throughout cabinet interior
- Available in three finish options (OSLO Cherry, OSLO Oak and OSLO Figured Mahogany)
- Specifications: 18" (D) x 21" (W) x 29.5" (H)

### PRODUCTS

### Description

NOBS3D-C	Oslo Bedside Cabinet, 3Drawer, Cherry
NOBS3D-FM	Oslo Bedside Cabinet, 3Drawer, FIG Mahogan
NOBS3D-O	Oslo Bedside Cabinet, 3Drawer, Oak