

Dual Oxygen Cart

HCPCS:




FIGURE A



- Lightweight silver vein cart can accommodate two "D" or "E" size cylinders
- Handle is adjustable to accommodate most individuals
- Silver vein finish easy-to-maintain & less institutional looking (Figure A)
- Large, composite wheels allow for easy, safe maneuvering

PRODUCTS

Description

| | |
|-----------|--------------------------------------|
| 13001SV-2 | Oxygen CylinderCart,Dual,Silver,2/cs |
|-----------|--------------------------------------|

SPECIFICATIONS

Description

| | |
|------------------------|------------------|
| Warranty | Limited Lifetime |
| Weight Capacity | 80 lbs |
| Height | 12.2" |
| Length | 38.6"-41.9" |
| Carton Shipping Weight | 22 lbs. |
| Weight | 5.4 lbs. |