

## Hamper Stand

HCPCS:



- Supports 36 42 gallon bags
- Glides easily with 3" locking casters
- Chrome-plated steel is attractive and easy-to-maintain
- Non-slip pedal raises the white, powder-coated steel lid
- Ships knocked down

### PRODUCTS

### Description

PRODUCTS	Description
13070	Hamper Stand, 1/cs

### SPECIFICATIONS

### Description

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Height	38"
Base	16.25" (W) x 19.25" (D)
Carton Shipping Weight	28 lbs.
Weight	24 lbs.