

## GRAVIS Floor Lift, Power Base 600 lb

HCPCS: E0635



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F



FIGURE G



FIGURE H

- Weight capacity of 600 lb
- Base option: Power
- Elegant, contemporary design with soft lines but still heavy duty
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Audible low battery alert fuel indicator, if maintenance is due, and if lift detects weight overload
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Maintenance monitoring: by months or lift cycles
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS	Description
FLP600	Gravis Floor lift, power base, 600 lbs
SPECIFICATIONS	Description
Battery Pack Weight	7.9 lb (5.0A), 4.9 lb (2.9A)
Boom Height	26" - 76"
Ground Clearance	2.9"
Warranty on Battery	1 yr
ACCESSORIES	Description
FLSCALE	Patient Lift Scale 660lbs capacity, EA
13012	Patient Sling,Canvas,1/bx, BX
13019	Chains for Patient Lft,1pr/bx, BX
13025	Sling,Polyst. for Lift 1/ea, BX
13026	Sling,Polyst. w/Comde Open1/ea, BX
13060	BariatricSling,Canvas,1/bx, BX
13061	BariatricSlingw/ComdeOpen,Canvas,1/bx, BX
13232	Straps for Patient Slings,1set, BX
13220L	U-Sling 4pt padded Lg.1/bx, BX
13220M	U Sling 4pt padded Med,1/bx, BX

**ACCESSORIES****Description**

13220S	U-Sling 4pt padded Sm,1/bx, BX
13220XL	U-Sling 4pt padded X-lg,1/bx, BX