

## Pivot & Tilt Overbed Table

HCPCS: E0274



FIGURE A

- The mast of the table pivots and can be locked in one of the three positions from flat on the floor to 90° (Figure A)
- The tabletop tilts with the angle of the mast to provide a flat surface regardless of the angle of the mast
- The pivot feature allows the table top to be positioned closer to individuals in a bed, sitting in a wheelchair, or folded for storage
- Four 2" casters (two locking)
- Only wheels have to be installed
- New and improved larger tubing

| PRODUCTS               | Description                   |
|------------------------|-------------------------------|
| 13000                  | Overbed Table,Pivot/Tilt,1/cs |
| SPECIFICATIONS         | Description                   |
| Warranty               | Limited Lifetime              |
| Weight Capacity        | 40 lb                         |
| Carton Shipping Weight | 22 lb                         |
| Tabletop Adjustment    | 28" - 37"                     |
| Weight                 | 21 lb                         |