

Deluxe Automatic Blood Pressure Monitor, Upper Arm

HCPCS: A4670-1,A4670-1



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- X-Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory stores 120 readings (2 x 60)
- Average Last 3 Results

PRODUCTS

Description

BP3400	BP Monitor, Auto, Upper Arm, 1/EA
--------	-----------------------------------

SPECIFICATIONS

Description

Warranty	2 Years
Measurement Range	Pressure: 0mmHg~300 mmHg Pressure: +/- 3 mmHg
Power Supply	4x AA Batteries or AC Adapter