

## Deluxe Automatic Blood Pressure Monitor, Wrist

HCPCS: A4670-1,A4670-1



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory Stores 120 Readings (2 x 60)
- Average Last 3 Results

### PRODUCTS

### Description

|        |                                      |
|--------|--------------------------------------|
| BP3200 | BP Monitor, Deluxe, Auto,Wrist, 1/EA |
|--------|--------------------------------------|

### SPECIFICATIONS

### Description

|                   |   |
|-------------------|---|
| Warranty          | 2 Years                                     |
| Measurement Range | Pressure: 0mmHg~300mmHg Pressure: +/- 3mmHg |
| Power Supply      | 2x AAA Batteries                            |