

PainAway Plus Wireless TENS

HCPCS:



FIGURE A



FIGURE B

- Small, portable and discreet
- Single-channel transcutaneous electrical nerve stimulator
- Can be used for chronic and acute pain on most areas of the body
- Wireless design with rechargeable receiver no batteries to purchase
- Long lasting electrode pad for up to 30 uses
- Provides constant current throughout each treatment 4 pulse modes automatically adjust
- Intensity Control is adjustable from 0-80 mA into 500 ohms 20 adjustable levels
- Amplitude range: maximum: 50V-90V minimum: 0-10V (1KO load)
- Frequency: 1.2-60Hz (±15%)
- Pulse: 100°S (±30%)
- Comes complete with Control Unit, 5"x3" Long Use Gel Electrode Pad, Recharging Kit and Instruction Manual
- Electrode refills available for separate purchase (RTLAGF-810 5"x3" electrode included)

PRODUCTS

Description

RTLAGF-800

PainAwayPlus Wireless TENS, 1EA/RTL

SPECIFICATIONS

Description

Warranty

1 Year on Receiver