

## Med-Aire Alternating Pressure Pump and Pad System

HCPCS:



- Quiet pump alternately inflates and deflates 130 individual bubble air cells
- 4 LPM (liters per minute) pump produces consistent air flow and pressure over a 5 minute cycle time and is available fixed or variable (Fig.A)
- Variable pressure setting on pump allows for individualized comfort increasing patient compliance
- Durable, heavy gauge, premium vinyl bubble pad provides superior therapeutic treatment and is available with or without end flaps
- Pump has built-in brackets for mounting easily to bed

PRODUCTS	Description
14001E	PP Sys,Dlx,Variable Pressure,1/ea 6/cv
14001EFD	Dlx APP w/Pad w/End Flaps,1/CV, 6/CT
14002E	APP Sys,Standard,Fixed Pressure,1/ea,RTL

SPECIFICATIONS	Description
Weight Capacity	300 lbs

ACCESSORIES	Description
14001T	Tubing for 14001/14002,1/ea, EA
14003	APP Pad Only, 1/ea, RTL, EA
14005E	Pump only 14001E, 1/EA, EA
14006E	Pump only for 14002E,1/ea, EA