

HCPCS:

| PRODUCTS   | Description                          |
|------------|--------------------------------------|
| STDS2A4300 | Rear wheel assembly,24",Cougar,1/ea  |
| STDS2A4301 | Armrest 'T' Frame,L3 Coug.1/ea       |
| STDS2A4302 | Arm Receiver,L3 Coug.1/ea            |
| STDS2A4304 | Caster Fork,L3 Coug.1/ea             |
| STDS2A4306 | Seat Uphol.16"L3Coug.1/ea            |
| STDS2A4307 | Back Cane,L3 Coug.1/ea               |
| STDS2A4308 | Seat Uphol.18"L3 Coug.1/ea           |
| STDS2A4309 | Footrest Rcvr.Lock,L3Coug.1/ea       |
| STDS2A4310 | Caster,8"L3 Coug.1/ea                |
| STDS2A4314 | Caster Bearing, 8mm x 22mm, 1/ea     |
| STDS2A4316 | Back Uphol.16"L3Coug.1/ea            |
| STDS2A4318 | Back Uphol.18"L3Coug.1/ea            |
| STDS2A431L | Arm Panel,Left,Cougar,1/ea           |
| STDS2A431R | Arm Panel,Right,Cougar,1/ea          |
| STDS2A4321 | Bearing Housing, Left, L3 Coug, 1/ea |
| STDS2A4324 | Quick release axle,Cougar,1/ea       |
| STDS2A432L | Brake, left, Cougar, 1/ea            |
| STDS2A432R | Brake, Right, Cougar, 1/ea           |
| STDS2A4333 | Fork Stem Bearing,L3Coug.1/ea        |
| STDS2A4335 | Brake Tip, L3 Coug.1/ea              |
| STDS2A4336 | Back Support Bar, Cougar, 1/ea       |
| STDS4058   | Seat Rail Saddle for Viper GT, 1/ea  |
| STDS871    | Hand Grip,Sentra,1/ea                |