

Super Absorbent Commode Pail Liners



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

• The new Drive Commode Liners are a single-use medical device to be positioned over a commode. These Commode Liners feature an absorbent powder which reduces splashing and eliminates odor. The new Drive liners are a convenient solution for cleaning and removing the waste from a commode. These highly absorbent liners can hold up to two quarts of liquid. The Commode Liners feature tie handles for easy and quick disposal. They come with a dispensing package for easy transport and storage. The new Commode Pail Liners provides comfort and dignity to the patient and patient's family!

- Holds up to 2 quarts or 64 ounces of liquids
- Avoid spills and splashes
- Commode liner can be used for multiple uses
- Easy to use, convenient and safe
- Equipped with a super absorbent powder which soaks up bodily waste (up to 64 oz) within seconds
- Prevents the soiling/contamination of commode pails and reduces unpleasant odors
- Provides comfort and dignity to the patient and their family
- Safe and secure collection, tie handle for quick and easy transport and disposal of contaminated bodily fluids
- Time and money saving: No hand washing
- Includes 7 liners per box, and 6 boxes per case
- Compatible with; RTL11148KDR, 11148N-4, 11148-4, 11148-1, 11148CE-4, 11105N-4, 11112-2, 11114KD-1, NRS185006, NRS185007, RTL11158KDR, 11125KD-1, 11125KD-2, 11125PSKD-1, 11120SV-1, 11125KD-1A, 11117N-2, 11117N-1, 11130-2, 11135-1, NRS185008

PRODUCTS

Description

RTL12095	Commode Liners, 7 Liners/BX, 6BX/CV
----------	-------------------------------------

SPECIFICATIONS

Description

Retail Box Dimensions	1.75" (W) x 7.5" (D) x 3.5" (Height without 5th panel) or 4.75" (Height with 5th panel)
-----------------------	---

Weight	0.69 lb
--------	---------