

P750 Trapeze



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- The free-standing bariatric trapeze, P750TPZ, is a floor-mounted device and optional accessory to the P750 LTC Bed. This heavy-duty trapeze is intended to support a portion of an individual's weight to aid in repositioning of their body while in bed, or assist in the transfer from bed to a chair.
- Increase resident independence with adjustable handle with strap and 5 positioning options (2" increments).
- Ensure long-term durability with a welded steel frame.
- Fits seamlessly with existing Drive DeVilbiss products using a wide, sturdy base designed for most beds.
- Saves time on room setup and storage with a wheeled base for easy transport.
- Weight Capacity: 600lb.
- Warranty: 2 Years, 10 Years Weld

PRODUCTS

Description

P750TPZ	Trapèze P750
---------	--------------

SPECIFICATIONS

Description

Garantie	2 Years, 10 Years Weld
Poids du produit	127lb.
Capacité de poids	600lb.
Dégagement de la base	1.5"
Longueur de la flèche	45"
Poids à l'expédition du carton	138lb.
Dimensions	75.25" (H) x 36.5 (W) x 59.25" (D)
Largeur entre les bras	29.25"