

Half-Length Bed Rail

HCPCS:




FIGURE A



FIGURE B




FIGURE C



FIGURE D

- Half rail was designed to comply as part of the FDA Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment
- Tool-free width adjustment allows the rail to be adjusted in 2" increments from 36" 42" (Figure A - B)
- Ideal for protection when bed is elevated
- Chrome finish on second and fifth bar allows for scratch-proof movement of adjustment bracket
- Fits most beds with steel frame (Figure C)
- Attractive brown-vein finish is easy to maintain and less institutional looking (Figure D)
- Constructed of 1" steel
- Not for use with "Pan-Style" bed surfaces
- Easy to install

PRODUCTS

Description

| | |
|---------|--|
| 15208BV | HalfRail,ToolFreeAdjWidth,BrownVein,1/cs |
|---------|--|

SPECIFICATIONS

Description

| | |
|------------------------|--------------------|
| Garantie | Limited Lifetime |
| Réglage de la longueur | 36", 38", 40", 42" |