

P Pod

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Provides midline alignment of trunk, pelvis, and lower extremities Simple/Lightweight design is easy for parents and caregivers to use and clean Adjustable angle based on alignment of Pod on bean bag base Provides clinically appropriate postural support through use of the molded foam liner Pelvic belt keeps user in desired position and headrest provides additional head positioning Available in 4 Sizes (Small/Medium/Large/ Extra-Large) and now available in wipe-clean Dartex fabric

PRODUCTS

Description

PRODUCTS	Description
APP-2000S-DX	Support P Pod, S, Dartex lavable
APP-2000S-MP	P Pod, S, Rose Mermaid
APP-2000S-PB	P Pod, S, bleu piscine
APP-3000M-DX	Support P Pod, M, Dartex lavable
APP-3000M-MP	Support P Pod, M, rose sirène
APP-3000M-PB	P Pod, M, bleu piscine
APP-4000L-DX	Support P Pod, L, Dartex lavable
APP-4000L-JP	P Pod, L, méduse violette
APP-4000L-MP	P Pod, L, Mermaid Rose
APP-4000L-WB	P Pod, L, bleu baleine
APP-5000XL-DX	Support P Pod, XL, Dartex lavable
APP-5000XL-JP	P Pod, XL, méduse violette
APP-5000XL-MP	P Pod, XL, Rose Mermaid
APP-5000XL-WB	P Pod, XL, bleu baleine