

## P Pod

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Provides midline alignment of trunk, pelvis, and lower extremities Simple/Lightweight design is easy for parents and caregivers to use and clean Adjustable angle based on alignment of Pod on bean bag base Provides clinically appropriate postural support through use of the molded foam liner Pelvic belt keeps user in desired position and headrest provides additional head positioning Available in 4 Sizes (Small/Medium/Large/ Extra-Large) and now available in wipe-clean Dartex fabric

| PRODUCTS      | Description                       |
|---------------|-----------------------------------|
| APP-2000S-DX  | Support P Pod, S, Dartex lavable  |
| APP-2000S-MP  | P Pod, S, Rose Meraid             |
| APP-2000S-PB  | P Pod, S, bleu piscine            |
| APP-3000M-DX  | Support P Pod, M, Dartex lavable  |
| APP-3000M-MP  | Support P Pod, M, rose sirène     |
| APP-3000M-PB  | P Pod, M, bleu piscine            |
| APP-4000L-DX  | Support P Pod, L, Dartex lavable  |
| APP-4000L-JP  | P Pod, L, méduse violette         |
| APP-4000L-MP  | P Pod, L, Meraid Rose             |
| APP-4000L-WB  | P Pod, L, bleu baleine            |
| APP-5000XL-DX | Support P Pod, XL, Dartex lavable |
| APP-5000XL-JP | P Pod, XL, méduse violette        |
| APP-5000XL-MP | P Pod, XL, Rose Meraid            |
| APP-5000XL-WB | P Pod, XL, bleu baleine           |