

Nitro Sprint Cane Holder

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

The Nitro Cane Holder is designed to fit on the Nitro Sprint Rollator with built-in accessory mount to provide convenience at every turn

- Fits on all Nitro Sprint Rollators with built-in accessory mount
- Securely holds most standard canes on side of rollator
- Easily mount cane or remove cane from rollator for seamless transition from cane to rollator

PRODUCTS

102662-CANE

Description

Porte-canne Nitro, 1BX

SPECIFICATIONS

dimensions du produit

Description

3,8" (L) x 2,4" (L) x 2" (H)

Prop. CA 65

Non