

## Nitro Sprint Cup Holder



HCPCS:



FIGURE A




FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F

The Nitro Sprint and Glide Cup Holder is designed to fit on all Nitro Sprint Rollators, Nitro Glide Knee Walkers, and MD by Drive Air-Go Walker with built-in accessory mount to provide convenience at every turn

- Fits on all Nitro Sprint Rollators, Nitro Glide Knee Walkers, and MD by Drive Air-Go Walker with built-in accessory mount
- Securely holds cups, water bottles, and travel mugs all at an arms reach

### PRODUCTS

### Description

102662-CUP	Nitro Universal Cup Holder,1EA
------------	--------------------------------

### SPECIFICATIONS

### Description

Product Dimensions	4.1" (L) x 4.3" (W) x 4.4" (H)
CA Prop 65	no