

ProBasics Round Handle Cane (Silver),

HCPCS:

- The ProBasics Round Handle Cane provides excellent support and is designed for patients who need additional assistance with balance.
- The round handle allows users to hook the cane over their arm while opening doors or performing other tasks.
- Supports patient weights up to 250 pounds.
- Easy one-button height adjustment with a locking anti-rattle ring allowing for quick, safe adjustment
- Adjusts in 1-inch increments
- Anodized aluminum finish with vinyl grip
- Curved handle allows cane to be retained on arm
- Easy-to-clean aluminum frame

PRODUCTS	Description
CNRSL	ROUND HANDLE CANE, SILVER 10EA/CV