

5L Compact Oxygen Concentrator

HCPCS:



- One of the smallest units on the market, the new 555 Oxygen Concentrator's compact size, durable construction, and quiet operation help you run more efficiently and reduce long-term ownership costs.
- Warehouse and fit more units per pallet or truckload with space-saving design
- Faster, safer deliveries at just 29.9 pounds, the 555 is 10 to 15 pounds lighter than competing concentrators
- Our most reliable concentrator ever with next-generation design leveraging decades of expertise
- Fewer potential opportunities for failure with reduced connections and clamps alongside simplified and streamlined assembly
- Quick servicing thanks to modular internal design providing easy access to streamlined components
- Extended component life due to reduced power consumption for lower internal temperatures
- Enhanced maneuverability and fewer tripping hazards with a compact, lightweight build
- Fewer noise complaints due to low and pleasant sound level
- Helps lower monthly energy costs through a 17% reduction in power consumption
- User peace of mind with easy-to-read digital controls that show unit status at a glance

PRODUCTS

Description

| | |
|-------|--------------------------------|
| 555DS | 5L Compact Oxygen Concentrator |
|-------|--------------------------------|
