



Padded Swivel Seat Cushion

HCPCS:



FIGURE A



FIGURE B

- Swivel seat allows individual to turn up to 360 degrees
- Makes getting on or off a seat or chair easy and comfortable (Figure A)
- Durable and stable
- Portable: weighs only 3.25 lbs
- Dimensions: 15 /" diameter
- Retail Packaged
- Weight Capacity 300 lbs

PRODUCTS

RTLAGEF-300

Description

Swivel Seat Cushion,Gray,RTL,1EA,8/CS