

## Pure Expression Carry Bag

HCPCS:



- Carry bag for RTLBP2000

| PRODUCTS | Description                             |
|----------|-----------------------------------------|
| BP001    | Breast Pump Carry Bag,1/ea              |
| BP0011   | PowerAdapter,6V,RTLBP2000/RTLBP2200, 1e |
| BP0012   | ExternalBatteryPack,6V,RTLBP2000/BP2200 |
| BP0013   | PumpKit w/o Insert,31MM,BP1000/2000/220 |
| BP0014   | Silicone Insert, 27MM, BP1000/2000/2200 |
| BP0015   | FeedingNippleSet, BP1000/2000/2200,1/pk |
| BP0016   | Base for Bottles, BP1000/2000/2200, 1/e |
| BP002    | Insulated Cooler Bag, 1/ea              |
| BP005    | PumpKit w/o Insert,27MM,BP1000/2000/220 |
| BP006    | Silicone Insert, 24MM, BP1000/2000/2200 |
| BP007    | ValveSetW/Membrane,BP1000/2000/2200,1/p |
| BP008    | Membranes ,BP1000/2000/2200, 4/bag      |
| BP009    | Tubing Set, BP1000/2000/2200, 1/pr      |