

CARRY CASE FOR SINGLE BED END

HCPCS:

- The new ProBasics Bed End Carry Case features strong wraparound handles and a durable polyester/PVC-coated shell with padded ends and base.
- It makes transporting 36" wide bed ends easier while protecting end caps from damage and panels from puncture or scratching during transit.
- Just bag the bed end as clean or dirty, slide the carry case over the covered bed end, close the hook-and-loop strap and you're ready to go!

PRODUCTS	Description
PBHBCASE	CARRY CASE FOR SINGLE BED END 25EA/CV