

AquaSense Adjustable Bath Seat with Backrest

HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Provide comfort and safety while showering or bathing. The AquaSense® Adjustable Bath Seats with Backrest are made from lightweight aluminum, with blow-molded plastic for the seat and the backrest. They feature a textured surface, large drainage holes, handles in both seat and back, and slip resistant rubber safety tips. Adjustable legs raise or lower the seat height.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 300 lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 14" to 21".
- Drainage holes help to prevent water from pooling, making the bench more comfortable, less slippery and ensuring transfer from bathing to exiting the tub is more secure.
- Textured seat surface helps to prevent slipping.
- Push button, snap together design requires no tools and assembles in just minutes.
- Feel good knowing that plastics parts of this product are made using recycled materials.

PRODUCTS

Description

770-510	AQS B.SEAT W/ BK WHITE, 1EA, 2CV, Ret.
---------	--

SPECIFICATIONS

Description

Weight Capacity	300 lb
-----------------	--------