

Shower Chair with Back and Removable Padded Arms

HCPCS:



FIGURE A



- Removable padded armrests provide comfort and added security
- Blow-molded bench and back provide comfort and strength
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Leg height adjusts in 1" increments
- Made using recycled materials

PRODUCTS	Description
12445KD-1	KD Bath Bench w/Back & Padded Arm,1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	300 lbs
Carton Shipping Weight	21 lb
Weight	9.75 lb
Width Between Arms	18.5"
ACCESSORIES	Description
12002T	Tips for 12002,4/pkg., BX