

## Niagara Bathing System

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- CONTACT YOUR SALES REP TO PLACE AN ORDER. This product is not purchasable online.
- Thermostatic mixing valve - Provides an anti-scald protection that dispenses 15 gallons per minute for fast tub fill
- Built-in disinfection system - The most time-efficient solution for disinfection of the baths external surfaces and minimizes the risk of cross-infection
- The built-in disinfection system cabinet is conveniently located on the side of the tub for easy access
- The lockable cabinet provides a secure place for cleaning supplies and disinfectants
- Handheld shower wand- Hair and upper body washing can be used at the same time the tub is filling
- Water temperature display - Monitors tub-fill and shower wand temperatures
- 3/4" supply lines - Constructed of stainless steel mesh clad
- The CleanDrive Air Spa System - Variable speeds and pulse modes designed to enhance the bathing experience for the resident and the caregiver

### PRODUCTS

### Description

DB-NIAGARA-1200	Niagara Bathing System Fixed Base
DB-NIAGARA-1300	Niagara Bathing System Fixed Air Spa
DB-NIAGARA-1500	Niagara Bathing System ADJ Base
DB-NIAGARA-1600	Niagara Bathing System ADJ Air Spa

### SPECIFICATIONS

### Description

Weight Capacity	600 lbs
Dimensions	Fixed height: 35.5" (W) x 84" (L) x 32.5" (H) Adjustable Height: 35.5" (W) x 90" (L) x 30"- 42" (H)