

## Saratoga Spring Bathing System

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Thermostatic mixing valve - Provides an anti-scald protection that dispenses 15 gallons per minute for fast tub fill
- Built-in disinfection system - The most time-efficient solution for disinfection of the baths external surfaces and minimizes the risk of cross-infection
- The built-in disinfection system cabinet is conveniently located on the back end of the tub for easy access
- The lockable cabinet provides a secure place for cleaning supplies and disinfectants
- Handheld shower wand - Hair and upper body washing can be used at the same time the tub is filling
- Water temperature display - Monitors tub-fill and shower wand temperatures
- 3/4" supply lines - Constructed of stainless steel mesh clad
- Three plumbing configurations - to choose from side, front, or island configuration
- The island configuration allows 360° unobstructed access for the caregiver and with no visible hook-ups or cords, the Saratoga Springs Bathing System will give the bathing area the look of a luxury spa
- The CleanDrive Air Spa System - Variable speeds and pulse modes designed to enhance the bathing experience for the resident and caregiver

PRODUCTS	Description
DB-SARATOGA-1000	Saratoga SE Bathing System Base
DB-SARATOGA-1100	Saratoga SE Bathing System Air Spa
DB-SARATOGA-1200	Saratoga SE Bathing System Auto Fill
DB-SARATOGA-1300	Saratoga SE Bathing System AS AF
DB-SARATOGA-SCALE	Saratoga Scale Add On

SPECIFICATIONS	Description
Seat Dimensions	24" (W) x 21" (H)
Weight Capacity	450 lb or Resident 730 lb or Resident and Water
Height	78"
Dimensions	34" (W) x 73" (L) x 51" (H)