

Carex White Bathtub Rail

HCPCS:



- The Carex White Bathtub Rail provides a firm and sturdy grip, allowing for safe entrance in and out of the tub.
- It features a wide contact area that allows even distribution of weight.
- Safely assists those who have difficulty getting in and out of the tub and prevents falls
- Fits most bathtubs (flat top walls 2.5 to 5 inches wide); ideal for use on fiberglass bathtubs
- Wide contact area for even weight distribution
- Attractive white finish
- Attaches easily with adjustable mounting clamp, no tools required
- Extends up to 14 inches above tub wall
- Weight capacity: 200 lbs
- Retail Packaging

PRODUCTS	Description
FGB20400 0000	WHITE BATHTUB RAIL 6EA/CV