

Carex E-Z Lock Raised Toilet Seat W/ Armrests

HCPCS:

- Our E-Z Lock Raised Toilet Seat with Armrests adds 5 inches of height to the toilet seat to assist those with bending or sitting difficulties.
- The easy-to-use locking mechanism secures the raised toilet seat to the toilet bowl and the smooth contour seat has built-in hand grips for easy handling.
- This seat fits most toilets, is lightweight, and easy to clean.
- Adds 5" of height to the toilet seat
- 17.25" of width between the armrests
- Comes with padded arm rests
- Secure locking mechanism included
- Retail Packaging

PRODUCTS	Description
FGB30400 0000	E-Z LOCK RAISED TOILT SEAT W/ARMS 2EA/CV