

## Carex Universal Transfer Bench, 2/CA

HCPCS:

- The Carex Universal Transfer Bench makes it easy for users to safely transfer into and out of the tub.
- It rests with one side of the bench in the tub and one side out, so users can simply sit on the outside and slide their legs around and into the bath.
- Exact level height adjustable legs (17.5" to 22.5") make it easy to level the bath bench to the perfect height
- Converts to either right- or left-hand entry
- Anti-slip feet provide safety and stability
- Tool-free assembly and disassembly for quick and easy set-up, storage and travel
- Features a backrest with two built-in handles to make positioning easier
- Fits most bathtubs and showers
- Made of sturdy, comfortable, and noncorrosive materials
- Weight capacity: 400 lbs
- Quantity: 2/case

PRODUCTS	Description
FGB17000 0000	UNIVERSAL TRANSFER BENCH 2EA/CV