



ProBasics Bariatric Commode with Extra Wide Seat

HCPCS:

- ProBasics Bariatric Commode Extra Wide Seat features a 24-inch wide platform seat and reinforced steel frame to support patient weights up to 650 pounds and is ideal for larger patients.
- It comes complete with commode bucket and lid.
- The versatile design includes a commode bucket with cover, armrests for added comfort, angled legs for stability and non-marring rubber tips
- Welded steel construction with 24-inch extra wide snap on platform seat designed for strength and stability
- Tool-free assembly
- Commode opening is contoured for increased comfort
- Seat height readily adjusts in one-inch increments
- Angled legs for stability with non-marring rubber tips
- Plastic armrests provide added comfort
- Comes complete with a commode bucket and cover
- Bucket capacity: 9.3 quarts / 8.8 liters

PRODUCTS

Description

BSB24C

BARIATRIC STEEL COMMUNE-24"W 2EA/CV