

## ProBasics Bariatric Shower Bench without Back

HCPCS:

- The ProBasics Bariatric Bath Bench is designed for individuals who require seating support and stability when bathing or showering.
- The bench features a special cross brace design with a heavy duty blow-molded textured seat that can support patient weights up to 550 pounds.
- The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety.
- Tool-free assembly and limited lifetime warranty are standard.
- Comfortable, blow-molded seat with integrated handles for additional support
- Dual cross brace design reinforced with aircraft-style rivets for extra stability
- High-strength aluminum frame to prevent corrosion
- Textured contoured surface with drainage holes to reduce slipping and increase comfort
- Wide-based tips help increase stability
- Legs are height adjustable in one-inch increments with heavy duty height adjustment pins to accommodate larger individuals
- Quick tool-free assembly

**PRODUCTS**

**Description**

---

BSBC	BARIATRIC SHOWER CHAIR W/O BACK 2EA/CV
------	--

---