

ProBasics Raised Toilet Seat with Lock

HCPCS:

- ProBasics Raised Toilet Seats are a must for anyone with bending or sitting difficulties.
- Constructed of heavy duty blow-molded plastic that is durable and easy to clean.
- Adjustable front locking mechanism tightens easily with one hand and fits most toilets.
- Prevents shifting of seat when transferring.
- Elevates existing seat height by 4.5"
- Extra-large hygienic hole cutout
- Tool-free assembly
- Color: white
- 350 lb weight capacity
- HCPCS Code: E0244*

PRODUCTS

Description

BSRTSL	RAISED TOILET ST W/LOCK 1EA/
--------	------------------------------
