

ProBasics Steel Drop-Arm Commode,

HCPCS:



- The ProBasics Steel Drop-Arm Commode features a durable steel frame with an easy-release push-button arm mechanism allowing for safe lateral patient transfers to and from the commode.
- The anti-slip, soft Hypalon® armrests provide extra comfort and support.
- The frame can be configured to be used as a bedside commode, toilet safety frame or raised toilet seat.
- The tool-free removable back makes configuring the commode easy for any application.
- Durable steel frame provides strength and stability
- Equipped with easy-to-release arm mechanism allowing for safe lateral transfers to and from the commode
- Push-button drop-arm allows easier side transfer
- Soft, slip-resistant Hypalon® covered armrests
- Seat height readily adjusts in one-inch increments
- Bucket slides easily into place
- Removable backrest

PRODUCTS

Description

BSDAC

DROP ARM COMMODE W/SPL GRD 1EA/