

ProBasics Toilet Safety Frame

HCPCS:



- The ProBasics Toilet Safety Frame provides a safe, stable aid for individuals who struggle with the up and down demands of toileting.
- The frame is designed to accommodate both standard and elevated toilet seats.
- It is easy to clean and the aluminum frame mounts securely onto the bowl with an adjustable bracket.
- The toilet safety frame can safely support patient weights up to 300 pounds.
- Frame is made of anodized, non-corrosive one-inch aluminum tubing
- Arms are height and width adjustable
- Versatile design with tool-free assembly quickly secures to most toilets
- Contoured waterfall armrests provide additional comfort and support while providing an excellent gripping surface
- Powder-coated aluminum bracket easily attaches frame to toilet
- HCPCS Code: E0243*
- ProBasics Toilet Safety Frame Informational Video-

PRODUCTS	Description
BSTF	TOILET SAFETY FRAME 4EA/CV