

AquaSense Raised Toilet Seat with Lid



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

Designed for individuals who have difficulty sitting down or standing up from the toilet so they can easily get on and off the toilet. The AquaSense® Raised toilet seat with lid makes it easier to get on and off the toilet by putting the toilet seat 4 inches higher. Ideal for anyone with mobility issues, this toilet seat is easy to clean and features quick and easy tool-free installation.

- Adding 4" inches of height to your toilet making it easier and safer to sit or stand
- Includes a lid for better hygiene and aesthetic
- Contoured to provide a safe and comfortable elevation
- Made from PE that is both lightweight and durable to support up to 300 lbs
- No Tools Required! Simply fasten the clamps onto the toilet bowl. When not in use easily remove the toilet seat by releasing the clamps
- Fits most standard toilets
- Feel good knowing that plastics parts of this product are made using recycled materials

PRODUCTS

Description

770-626	AQS 4" RTS, LID, WHITE, 1EA,4CV,Ret
770-629	AQS Elon 4" RTS,LID, WHITE, 1EA,4CV,Re

SPECIFICATIONS

Description

Weight Capacity	300 lb (136 Kg)
-----------------	-----------------