

Premium Plastic Raised, Regular/Elongated Toilet Seat, with Lock

HCPCS: E0244



FIGURE A



FIGURE B

- Designed for individuals who have difficulty sitting down or standing up from the toilet
- Wide opening in front and back for personal hygiene
- Lightweight and portable
- Standard locking mechanism ensures safety (Figure A)
- Fits most elongated toilets
- No tools required for installation
- Easy-to-clean
- Made using recycled materials

| PRODUCTS | Description |
|------------------------|--------------------------------------|
| 12014 | Raised Toilet St. w/o Arms,1/cs, RTL |
| SPECIFICATIONS | Description |
| Warranty | Limited Lifetime |
| Seat Dimensions | 17" (W) x 16.5" (D) x 4.5" (H) |
| Weight Capacity | 300 lb |
| Carton Shipping Weight | 6 lb |
| Weight | 4 lb |