

PreserveTech Secure Lock Raised Toilet Seat

HCPCS: E0244,E0244



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F



FIGURE G

- The new PreserveTech Secure Lock Raised Toilet Seat features three locking points on the front and sides for added security and stability. This toilet seat adds an extra 5 of height to the toilet and includes padded armrests that can be adjusted or removed. Most non-permanent raised toilet seats only offer one or two locking points of security resulting in instability when the end user places more weight on one side of the unit. The seat fits on most elongated and standard toilets and installs without tools.
- Three locking points of security located on the front, left and right side providing the highest level of stability and security. Most raised toilet seats only offer one or two locking points of security.
- Adding 5 inches of height to your toilet making it easier and safer to sit or stand.
- Padded armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing up.
- Easily adjust armrest height to your needs.
- Contoured to provide a safe and comfortable elevation.
- Made from HDPE that is both lightweight and durable to support up to 350 lbs.
- Fits comfortably on most standard and elongated toilets.
- No Tools Required! Simply turn the adjustable knobs to tighten the locking plate and ensure a secure fit.
- Made using recycled materials.
- Drive PreserveTech products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.
- *Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

PRODUCTS

Description

RTL12C003-WH	PreserveTech Secure Lock RTS, 1EA, 1/CV
--------------	---

SPECIFICATIONS

Description

Weight Capacity	350 lb
Product Type	Toilet Seat
Width Between Arms	17.5"