

## Aquasense Adjustable Transfer Bench

HCPCS: E0247,E0247



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Slide effortlessly over the tub wall of your bathtub with your AquaSense® Transfer Bench.
- The Adjustable Bathtub Transfer Bench, by AquaSense®, is knocked-down and blow molded. It is built with a light weight anodized rust-resistant frame. The transfer bench's legs have nine (9) adjustments in 1.2" increments to accommodate different bath tub heights.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 250 lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- A removable push-button armrest supports transferring in and out of the bathtub.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 18" to 22".
- Drainage holes help to prevent water from pooling, making the bench more comfortable, less slippery and ensuring transfer from bathing to exiting the tub is more secure.
- Push button, snap together design requires no tools and assembles in just minutes.
- Feel good knowing that plastics parts of this product are made using recycled materials.

### PRODUCTS

### Description

770-406

AQUASENSE TRANSFER BENCH, 1EA, 2CV, Ret

### SPECIFICATIONS

### Description

Weight Capacity

250 lb (113 kg)

Seat To Floor Height

18 to 22 (45.7 to 55.9 cm)

Height Adjustment

4 in 0.47 (1.2 cm) increments

Overall Width

26 (66 cm)

Retail Box Dimensions

20.9 (L )x 17.7 (W) x 19.9 (H) |53 cm x 45 cm x 50.5 cm