

Bariatric Transfer Bench

HCPCS: E0248



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly of back, legs and arm
- A-frame provides additional stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- All transfer benches are reversible to accommodate any bathroom
- Euro-style clip ensures a secure fit (Figure A)
- No exposed hardware to injure patient
- Large, locking suction cups for added safety (Figure B)
- Weight Capacity 500 lbs
- Made using recycled materials

PRODUCTS

Description

| | |
|-----------|--------------------------------------|
| 12025KD-1 | Transfer Bench, Heavy Duty, KD, 1/cs |
|-----------|--------------------------------------|

SPECIFICATIONS

Description

| | |
|------------------------|-----------------------|
| Warranty | Limited Lifetime |
| Seat Dimensions | 26.5" (W) x 16.5" (D) |
| Weight Capacity | 500 lb |
| Carton Shipping Weight | 27 lb |
| Weight | 14 lb |