

Folding Universal Sliding Transfer Bench

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Transfer Bench slides left and right along frame for easy entry into and out of tub
- Seat height adjusts in 1/2" increments
- Folds flat for transport and storage (Figure A)
- Comes standard with removable soap dish (Figure B)
- Comes with safety net to catch dropped personal items and prevent bending
- Comes with seat belt for added safety
- Retail packaged (Figure D)
- Made using recycled materials

PRODUCTS

Description

RTL12075	Sliding Shower Chair w/Soap Box,1EA,2/CS
----------	--

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	22.5" (W) x 18.5" (D)
Weight Capacity	300 lb
Height	21" - 25"
Carton Shipping Weight	38 lb
Weight	16 lb

ACCESSORIES

Description

12075-A	Arm, w/hardware, RTL12075,1ea, EA
---------	-----------------------------------